

# Medical Illumination System Two LED Trio | XLDS-S23



*Stock Photo - Options May Vary*

## FACTORY NEW

Medical Illumination System Two LED Trio: XLDS-S23. The energy-efficient System Two LED consumes half the electricity of its halogen counterpart. In addition, long-life 50,000-rated LEDs reduce the need to purchase and store replacement bulbs and eliminate procedure room downtime caused by unexpected bulb burn-outs.

Moreover, the output of the XLDS-S23 is rated at 120,000 Lux. Additionally, the XLDS-S23 features a Kelvin color temperature of 4300°, superior shadow reduction, and an unmatched MI arm system for optimal surgical site visibility and positioning.

System Two LED technology reduces heat, improving comfort for patients and staff compared to traditional halogen lights.

## FEATURES

- Green Technology
- Consumes half the electricity of its halogen counterpart
- Long-life 50,000-rated LED's
- Output rated at 130,000 Lux with a Kelvin color temperature of 4300°
- Superior shadow reduction and designed

## SPECIFICATIONS

- Light Head Dimensions: 25.2" Diameter x 4.75" Deep
- Arm (Articulated): 4.0" Diameter (tapered) x 26.0" Long
- Ceiling Casting: 17.0" Diameter x 4.7" Deep
- Ceiling Rod: 2.5" Diameter x 9.0"-43.0" Long (depending on ceiling height)
- Trio Ceiling Mount: Approx. 155.5 lbs

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)