

Medical Illumination MI-1000 LED Surgical Light Single Ceiling | XLD-SC



Stock Photo - Options May Vary

FACTORY NEW

Medical Illumination MI-1000 LED Surgical Light Single Ceiling Mount. The MI-1000 LED Surgical Light from Medical Illumination features exceptional quality, energy-efficient operation, and cool, bright output.

It provides a white, bright light from multiple LED spots, thus enabling excellent shadow control. Additional Features: The MI-1000 has an improved level of tissue rendition for your exams or minor procedures.

Furthermore, our lighting equipment (excluding replacement parts) has a three-year limited warranty. Long-life LEDs rated up to 50,000 hours are convenient as there are no replacement bulbs or fuses necessary, resulting in lower overall maintenance costs.

FEATURES

- Very Low Heat
- High Efficiency
- Half the Electric Consumption of Standard Halogen Lights
- Long Life L.E.D.s
- 4300 Kelvin Color Temperature

SPECIFICATIONS

- Light Intensity at 1 MTR: 100,000 lux
- 4300° Kelvin Color Temperature
- Rated Bulb Life: 50,000 hours
- Field Size at 1 MTR: 6" Diameter
- Diameter: 18" diameter, Multi-Faceted

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)