

Medical Illumination MI-1000 LED Surgical Light Portable Unit | XLD-FM



Stock Photo - Options May Vary

FACTORY NEW

The Medical Illumination MI-1000 LED Portable Surgical Light features exceptional quality, energy-efficient operation, and cool, bright output. The MI-1000 provides a white, bright light from multiple LED spots, thus enabling excellent shadow control and an improved level of tissue rendition for your exams or minor procedures.

Additionally, the MI-1000 features reliable Lighting equipment and comes with a 3-year MFG limited warranty. Furthermore, the MI-1000 includes Long-life LEDs that last up to 50,000 hours, eliminating bulb replacements and reducing maintenance costs. With versatile mounting and advanced technology, this portable surgical light delivers excellent value in modern medical lighting.

Moreover, for reference, the model number is XLD-FM.

FEATURES

- Very Low Heat
- High Efficiency
- Half the Electric Consumption of Standard Halogen Lights
- Long Life L.E.D.s
- 4300 Kelvin Color Temperature

SPECIFICATIONS

- Light Intensity at 1 MTR: 65,000 lux
- >4300° Kelvin Color Temperature
- Rated Bulb Life: 50,000 hours
- Field Size at 1 MTR: 6" Diameter
- Diameter: 18" diameter, Multi-Faceted

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)