

Medical Illumination MI-750 LED Surgical Lights Dual Ceiling | XLDP-DC



Stock Photo - Options May Vary

FACTORY NEW

Medical Illumination MI-750 LED Surgical Lights Dual Ceiling: XLDP-DC. The MI-750 Series is the perfect light for minor procedures and examinations in smaller spaces.

Furthermore, being our newest minor surgery/procedure light, the MI-750, features exceptional quality, energy-efficient operation, and cool, bright output.

MI-750 LED Dual Ceiling Model XLDP-DC Features: The MI-750 provides a white, bright light from multiple LED spots, thus enabling excellent shadow control and an improved level of tissue rendition for your exams or minor procedures. Furthermore, the MI-750 offers versatile mounting and outstanding value for advanced medical lighting technology.

Moreover, the MI-750 Series is available in a variety of models.

FEATURES

- Dual Ceiling Model: MI750-DC
- Very low heat
- Output rated at 50,000 LUX at 39.4" (1 meter) distance
- 4,300 degrees Kelvin color temperature
- High efficiency with half the electric consumption of standard halogen lights

SPECIFICATIONS

- Note: For complete Specifications click below to view.
- Specifications for MI 550/750/1000 (click to view)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)