

Medical Illumination MI-750 LED Surgical Lights Portable Mount | XLDP- FM



Stock Photo - Options May Vary

FACTORY NEW

Medical Illumination MI-750 LED Surgical Lights Portable Mount: XLDP-FM. The MI-750 Series is the perfect light for minor procedures and examinations in smaller spaces. Our newest minor surgery/procedure light, the MI-750 mobile surgical light, features exceptional quality, energy-efficient operation, and cool, bright output.

Additionally, the MI-750 LED Portable Model features bright white light from multiple LED spots. In addition, this design enables excellent shadow control and improves tissue rendition for exams or minor procedures. Furthermore, the XLDP-FM is unquestionably the best value for those seeking the latest technology in medical lighting.

Moreover, the MI-750 Series is available in a variety of models. Have Questions about the XLDP-FM?

FEATURES

- Portable/Floor Model: MI750 Portable / Mobile Unit
- Very low heat
- Output rated at 75,000 LUX at 39.4" (1 meter) distance
- 4,300 degrees Kelvin color temperature

• Color Rendering Index of 97

SPECIFICATIONS

- Note: For complete Specifications click [below](#) to view.
- Specifications for MI 550/750/1000 (click to view)

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com