

Medical Illumination MI-750 LED Surgical Lights Wall Mounted | XLDP- WM



Stock Photo - Options May Vary

FACTORY NEW

Medical Illumination MI-750 LED Surgical Lights Wall Mounted: XLDP-WM. The MI-750 wall-mounted procedure light is the perfect light for minor procedures and examinations in smaller spaces. Furthermore, being our newest minor surgery/procedure light, the MI-750 features exceptional quality, energy-efficient operation, and cool, bright output.

Additionally, the MI-750 LED Wall Mount Model features bright white light from multiple LED spots. This design enables excellent shadow control and improves tissue rendition for exams or minor procedures. With a variety of mounting options, the MI-750 is unquestionably the best value available for those seeking the latest technology in medical lighting.

Moreover, the MI-750 Series is available in a variety of models.

FEATURES

- Wall Mount Model: MI750-WM
- Very low heat
- Offers 75,000 LUX of light intensity
- 4,300 degrees Kelvin color temperature
- High efficiency with half the electric consumption of standard halogen lights

SPECIFICATIONS

- Note: For complete Specifications click below to view.
- Specifications for MI 550/750/1000 (click to view)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)