

Medical Illumination MI-550 LED Surgical Lights Dual Ceiling | XLDE-DC



Stock Photo - Options May Vary

FACTORY NEW

The MI-550 LED Surgical Light delivers optimal illumination for surgical procedures, examinations, and clinical work. It produces 55,000 LUX at 24 inches with 4,300°K pure white light and a Color Rendering Index of 97, providing vivid and accurate color visualization with exceptional shadow control. The LED technology generates virtually no heat, creating a comfortable working environment while consuming only half the power of standard halogen lights.

The compact light head measures 9.2 inches in diameter and weighs just 5.6 pounds, fitting seamlessly into existing surgical suites. The articulating arm offers extensive positioning flexibility with $\pm 15^\circ$ vertical movement, approximately 540° horizontal rotation, and multiple articulation points for precise light placement. The durable LED design requires minimal maintenance compared to traditional lighting systems.

FEATURES

- Output rated at 55,000 LUX at 24" (610mm) distance
- 4,300° K pure white illumination
- CRI (Color Rendering Index) of 97
- Near-zero heat from both light beam and light head housing
- High efficiency with half the electric consumption of standard halogen lights

SPECIFICATIONS

- Output: 55,000 Lux @ 24"
- Light Head Power Consumption: 16W
- Light Head Dimensions: 9.2" Dia x 2.6" H
- Supply Circuit: 100-240 VAC, 50/60Hz
- # of LED's: 8

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)