

# Burton AIM LED Double Ceiling | ALEDDC



Stock Photo - Options May Vary

## FACTORY NEW

The Burton ALEDDC AIM LED Double Ceiling is a diagnostic examination light featuring state-of-the-art LED lighting technology. It delivers 45,000 Lux at 1 meter with a 92 CRI rating for accurate color visibility essential to precise diagnosis. The light provides 10% more illumination than the AIM-50 model while using 55% less energy.

The AIM LED features a 20,000-hour lamp life, approximately ten times longer than traditional lights, significantly reducing bulb replacement costs. The unique Y-shaped design facilitates optimal positioning around the physician's head, with a removable sterilizable handle for adjusting light patterns and multiple tilting reflectors.

The ceiling-mounted arm system offers 360° limitless rotation with a maximum reach of 63 inches and a 10-foot diameter coverage area. The light provides adjustable focus and a variable light field diameter of 13 to 15 feet, with 54 inches of depth illumination for comprehensive surgical and diagnostic applications.

## FEATURES

- Efficient LED Light Source
- On-Hand Quick Release for Removable Handle
- Removable Sterilizable Handle to Adjust Light Pattern
- Larger and Adjustable Light Pattern
- Unique Y-shaped design facilitates placement of lighthouse around the Physician's head
- 360° Limitless arm-system-rotation on single Ceiling Models
- The Arm System provides a maximum reach of 63" and a diameter of 10 feet
- 5-year limited warranty

## SPECIFICATIONS

- Illuminance: 45,000 Lux at 1 Meter
- Color Temperature: 4,300°K
- CRI (Color Rendering Index): 92
- The diameter of the Lighthouse: is 20"
- Rated Life Of LED Lamp: 20,000 Hours
- Swivel Radius of Lamp Housing-Ceiling Mounted: 63"
- Power: 84 Watts
- Minimum Ceiling Height/Maximum Ceiling Height: 8ft/12ft

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)