

Medical Illumination System Two LED Trio w/HD Video Camera | XLDS-S23VC



Stock Photo - Options May Vary

FACTORY NEW

Medical Illumination System Two LED Trio w/HD Video Camera: XLDS-S23VC. GREEN TECHNOLOGY: The energy-efficient System Two LED consumes half the electricity of its halogen counterpart.

Furthermore, the XLDS-S23VC has Long-life 50,000-rated LEDs alleviate the need to purchase and store replacement bulbs and eliminate procedure room downtime due to unexpected bulb burn-outs. Additionally, the XLDS-S23VC features a 120,000 Lux output at a 4300° Kelvin color temperature for superior visibility and shadow reduction.

Additionally, its MI arm system ensures precise positioning. The LED cool technology reduces heat, enhancing comfort for both patients and staff.

FEATURES

- Green Technology
- Consumes half the electricity of its halogen counterpart
- Long-life 50,000-hour rated LEDs
- Output rated at 130,000 Lux with a Kelvin color temperature of 4300°
- Superior shadow reduction and designed

SPECIFICATIONS

- Light Head Dimensions: 25.2" Diameter x 4.75" Deep
- Arm (Articulated): 4.0" Diameter (tapered) x 26.0" Long
- Ceiling Casting: 17.0" Diameter x 4.7" Deep
- Ceiling Rod: 2.5" Diameter x 9.0"-43.0" Long (depending on ceiling height)
- Trio Ceiling Mount Weight: Approx. 155.5 lbs

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)