

# Medical Illumination System Two LED Duo w/HD Video Camera | XLDS-S2VC



Stock Photo - Options May Vary

## FACTORY NEW

Medical Illumination System Two LED Duo w/HD Video Camera: XLDS-S2VC. This LED surgical light delivers up to 130,000 lux for clear, bright viewing. In addition, it reduces shadows and spreads light evenly across the field. Surgeons see details with ease during every procedure.

The built-in HD camera also supports video recording and training. The XLDS-S2VC surgical light helps lower energy use and cuts costs over time. It uses less power than older halogen lights while still producing strong light.

Moreover, the long-lasting LEDs reduce the need for frequent replacements. The adjustable arm allows easy positioning for better control.

## FEATURES

- Green Technology
- Consumes half the electricity of its halogen counterpart
- Long-life 50,000-hour rated LEDs
- Output rated at 130,000 Lux with a Kelvin color temperature of 4300°
- Superior shadow reduction and designed

## SPECIFICATIONS

- Light Head Dimensions: 25.2" Diameter x 4.75" Deep
- Arm (Articulated): 4.0" Diameter (tapered) x 26.0" Long
- Ceiling Casting: 17.0" Diameter x 4.7" Deep
- Ceiling Rod: 2.5" Diameter x 9.0"-43.0" Long (depending on ceiling height)
- Solo Light Head Weight: Approx. 73.5 lbs

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)