

# Luvis S200 LED Minor Surgical Light | Single Ceiling Mounted



*Stock Photo - Options May Vary*

## **FACTORY NEW**

The Luvis S200 Single Ceiling LED Minor Surgical Light is perfect for minor procedures and examinations in small spaces. Additionally, the S200 provides a shadow-less effect that offers less shadow, even under a 50% screened condition.

Moreover, the fanless design generates no heat or noise. Furthermore, the S200 LED Operating Light features an adjustable size pattern of 18cm to 30cm.

The Luvis S200 LED is also available in a dual-head configuration or on a mobile stand.

## **FEATURES**

- Excellent Illumination
- Ergonomic & Friendly Design
- Detachable & Autoclavable Main Handle
- No UV Emission & Cool Light
- Excellent Depth of Illumination

## **SPECIFICATIONS**

- LEDs: 18
- CRI: 90
- Color Temperature: 3,800K /4,300K/4,800K
- Head Size: 13.78" x 13.39"
- Focus Size: 7" (Circle)/11.8" (Circle)

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)