

Luvis S300 LED Surgical Light | Single Ceiling Mounted



Stock Photo - Options May Vary

FACTORY NEW

The Luvis S300 Single Ceiling Mounted LED Surgical Light delivers shadowless illumination for surgical procedures, examinations, and dental applications. It provides up to 130,000 lux of light output with three selectable operational modes, each offering unique intensity and light patterns optimized for specific clinical needs.

The S300 features adjustable light field sizes, 3-step color temperature selection (3,800K, 4,300K, 4,800K), and a high Color Rendering Index of 96 for accurate tissue visualization. The detachable and autoclavable sterile handle ensures infection control compliance in sterile environments.

The system incorporates 42 LEDs with a 60,000-hour lifespan and operates on standard 100-240V AC power. Its lightweight suspension system and various mounting options—single ceiling, dual ceiling, wall mount, or mobile stand—provide flexibility for different surgical suite configurations.

FEATURES

- 3 Selectable Modes (Surgical/Dental/Exam)
- Detachable & Autoclavable Main Handle
- Intensity: 130,000 lux
- Excellent Scalytic Effect
- 3 Selectable Color Temperature Options

SPECIFICATIONS

- Intensity: Max. 130,000 lux (Surgery Setting)
- Color Temperature: 3,800K /4,300K/4,800K
- CRI: 96
- LED Life Span: 60,000 hours
- Input Power: 100-240Vac, 50/60Hz

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)