

Luvis M200 LED Surgical Light | Single Ceiling Mounted



Stock Photo - Options May Vary

FACTORY NEW

The Luvis M200 Single Ceiling LED Surgical Light System is perfect for any procedure. Furthermore, the M200 features exceptional quality with an energy-efficient operation.

In addition, some great features of the Luvis M200 LED Surgical Light are: Excellent Illumination Ergonomic & Friendly Design Detachable & Autoclavable Main Handle No UV Emission & Cool Light Excellent Depth of Illumination 4,300k Color Temperature Additionally, the 54 LEDs in each light head produce up to 120,000 lux in a pattern that can be adjusted from 7.8" to 11.8". The Luvis Surgical Light includes a synchronized control panel on the sub-handle for easy access during surgery.

FEATURES

- Excellent Illumination
- Ergonomic & Friendly Design
- Detachable & Autoclavable Main Handle
- No UV Emission & Cool Light
- Excellent Depth of Illumination

SPECIFICATIONS

- LEDs: 54
- CRI: 95
- Color Temperature: 3,800k/4,300k/4,800k
- Head Size: 19.6"
- Focus Size: 7.8" (Circle)/11.8" (Circle)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)