

Luvis M300 LED Surgical Light | Single Ceiling Mounted



Stock Photo - Options May Vary

FACTORY NEW

The Luvis M300 Single Ceiling Mounted LED Surgical Light System is perfect for any procedure. The M300 features exceptional quality with an energy-efficient operation. Furthermore, the 84 LEDs in the Lighthouse produce up to 160,000 lux.

The light pattern can be adjusted from 7.8" to 11.8". Additionally, the Luvis Surgical Light also features: Head-mounted control panel Cool light with no UV emissions High Color Rendering Index of 95 Smart Sensors automatically detect obstacles in the light field and increase the intensity to compensate for the reduction in light. The Luvis M300 Single Ceiling LED Surgical Light comes in single, dual, or mobile configurations.

It also includes a 2-year manufacturer's warranty.

FEATURES

- Excellent Illumination
- Excellent Dilution Effect
- Ergonomic & Friendly Design
- Detachable & Autoclavable Main Handle
- No UV Emission & Cool Light

SPECIFICATIONS

- LED's: 84
- CRI: 95
- Color Temperature: 3,800k/4,300k/4,800k
- Head Size: 23.6"
- Focus Size: 7.8" (Circle)/11.8" (Circle)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)