

Medical Illumination MI-550 Series LED Single Ceiling | XLDE-SC



Stock Photo - Options May Vary

FACTORY NEW

Medical Illumination XLDE-SC MI-550 Series LED Single Ceiling. Whether you're performing surgery, procedures, or exams, the lighting conditions must be optimal to visualize your work accurately. Nothing less is acceptable.

At Medical Illumination, our sole focus is making that happen. The MI Series Features: Lights that are significantly brighter, whiter, and cooler than any previous technology.

So you get the most vivid, accurate colors with remarkable shadow control — and virtually no heat emission — which makes our lights more comfortable to work under than any other lights. Moreover, the compact but powerful MI-550 is the perfect fit and offers the greatest value in an examination light.

FEATURES

- Output rated at 55,000 LUX at 24" (610mm) distance
- 4,300° K pure white illumination
- CRI (Color Rendering Index) of 97
- Whiter, more color-correct light than halogen exam lights
- Near-zero heat from both light beam and light head housing

SPECIFICATIONS

- Number of LED's: 8
- Output: 55,000 Lux @ 24"
- Focal Length: 24"
- Dimming Levels: 2
- Dimming Percentages: 100%, 50%

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)