

Medical Illumination System Two MH | S2SMH



Stock Photo - Options May Vary

FACTORY NEW

Medical Illumination System Two High-Intensity Discharge Surgery Light MH Series: S2SMH. The System Two MH series represents a new, unparalleled level of surgery lighting technology, known as High-Intensity Discharge (HID).

Compared to traditional halogen surgery lights, the MH series provides cool, white light with virtually zero heat. System Two MH Features: This all-new design heralds a notable increase in efficiency over traditional halogen light sources, producing over 100,000 lux output with only a 39-watt H.I.D.

bulb, versus a 150-watt halogen bulb! The heat has been virtually eliminated, while color temperature has been significantly increased to over 4,300° Kelvin.

FEATURES

- 102,000 lux (9,500 ft candles) output at 1 meter distance
- 4,300° Kelvin color temperature
- Fully adjustable light pattern (from 6.5" to 9")
- Durable ceiling casting for optimum stability, balance, and function
- 360° rotational contacts at all key points

SPECIFICATIONS

- Weight, Solo ceiling mount assembly: Approximately 73.5 lbs
- Dimensions: Ceiling Casting - 17.0" diameter x 4.7" deep
- Ceiling Rod: 2.5" diameter x 9.0" - 43.0" long (depending on ceiling height)
- Arm (articulated): 3.5" diameter (tapered) x 26.0" long
- Light Head Assembly: 23.0" diameter x 10.0" deep

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)