

Luvis S200 LED Minor Surgical Light | Portable Mount



Stock Photo - Options May Vary

FACTORY NEW

The Luvis S200 Portable LED Minor Surgical Light delivers 60,000 lux of high-intensity illumination. Consequently, clinicians achieve exceptional visibility for minor surgeries and exams. The 95 CRI rating ensures accurate color rendition for precise tissue visualization.

Efficiency powers this advanced LED light. For instance, the Luvis S200 consumes minimal energy compared to standard halogen units. Because the LEDs last 50,000 hours, you avoid the cost of frequent bulb replacements.

This durability maximizes your equipment uptime. Furthermore, the cool LED technology emits no UV radiation, which keeps the surgical field comfortable for staff. Control stays effortless and precise.

FEATURES

- Excellent Illumination
- Ergonomic & Friendly Design
- Detachable & Autoclavable Main Handle
- No UV Emission & Cool Light
- Excellent Depth of Illumination

SPECIFICATIONS

- LED's: 18
- CRI: 90
- Color Temperature: 3,800K /4,300K/4,800K
- Head Size: 13.78" x 13.39"
- Focus Size: 7" (Circle)/11.8" (Circle)

**Brochure and images are a general overview of the product. Actual product options may vary.*