

# Luvis S300 LED Surgical Light | Portable Mount



Stock Photo - Options May Vary

## FACTORY NEW

The Luvis S300 Portable LED Surgical Lights are perfect for surgery, procedures, and examinations. The S300 portable operating light system provides a shadowless effect that offers less shadow, even under a 50% screened condition.

The S300 boasts a maximum of 130,000 lux of light output, with three selectable modes for surgery, exam, and dental. Each mode is unique in the pattern and intensity displayed.

The S300 LED Operating Light also features: Adjustable size pattern (Surgical/Dental/Exam) 3-step color temperature selection Detachable and autoclavable sterile handle See Gallery Image to illustrate the on-board control panel and the features offered. The Luvis S300 LED is also available in a dual-head or single-head configuration.

## FEATURES

- Excellent Illumination
- Easily Transported on a Wheeled Base
- Contemporary Design
- Ergonomic & User-Friendly Design
- Detachable & Autoclavable Main Handle

## SPECIFICATIONS

- LEDs: 42
- CRI: 96
- Color Temperature: 3,800K /4,300K/4,800K
- Head Size: 17" Head
- Focus Size: 7" (Circle)/11.8" (Circle)

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)