

Luvis M200 LED Surgical Light | Portable Mount



Stock Photo - Options May Vary

FACTORY NEW

The Luvis M200 Portable LED Surgical Light System is perfect for any procedure.

The M200 features exceptional quality with an energy-efficient operation.

Some great features of the Luvis M200 LED Surgical Light are: Excellent Illumination Ergonomic & Friendly Design Detachable & Autoclavable Main Handle No UV Emission & Cool Light Excellent Depth of Illumination 4,300k Color Temperature The 54 LEDs in each M200 light head produce up to 120,000 lux in a pattern that can be adjusted from 7.8" to 11.8". The Luvis M200 Surgical Light also features a control panel synchronized with the sub-handle, allowing the user to easily access it during surgery.

FEATURES

- Excellent Illumination
- Ergonomic & Friendly Design
- Detachable & Autoclavable Main Handle
- No UV Emission & Cool Light
- Excellent Depth of Illumination

SPECIFICATIONS

- LEDs: 54
- CRI: 95
- Color Temperature: 3,800k/4,300k/4,800k
- Head Size: 19.6"
- Focus Size: 7.8" (Circle)/11.8" (Circle)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)