

Burton AIM LED Floorstand w/Casters | ALEDFL



Stock Photo - Options May Vary

FACTORY NEW

Burton ALEDFL AIM LED Floorstand w/Casters. Introducing Burton's NEW AIM LED Model ALEDFL, the most exciting revolution in diagnostic examination lighting. The AIM LED uses the latest advances in state-of-the-art lighting technology to achieve incredible performance, unsurpassed value, and energy efficiency.

Additional AIM LED Features: The AIM LED provides 10% more illumination (45,000 Lux at 1 meter) than even the high-performing AIM-50®. A 92 CRI rating produces very accurate color visibility, which is essential for a precise diagnosis.

With ten times greater (20,000 hours) lamp life than traditional light source products, you will save more than \$1,000 in bulb replacements. Additionally, the AIM LED uses 55% less energy.

FEATURES

- Efficient LED Light Source
- On-Hand Quick Release for Removable Handle
- Removable Sterilizable Handle to Adjust Light Pattern
- Stand-By On/Off Switch
- Larger and Adjustable Light Pattern

SPECIFICATIONS

- Illuminance: 45,000 Lux at 1 Meter
- Color Temperature: 4,300°K
- CRI (Color Rendering Index): 92
- The diameter of the Lighthead: is 20"
- Light Field Diameter: 13' - 15"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)