

Medical Illumination MI-550 Series LED Portable Mount | XLDE-FM



Stock Photo - Options May Vary

FACTORY NEW

Medical Illumination MI-550 Series LED Portable Mount. Whether you're performing surgery, procedures, or exams, the lighting conditions must be optimal so you can visualize your work with complete accuracy. Nothing less is acceptable.

At Medical Illumination, our sole focus is making that happen. The MI Series Features: Lights that are significantly brighter, whiter, and cooler than any previous technology.

Get vivid, accurate colors, excellent shadow control, and virtually no heat, making our lights more comfortable than others. All at a price you won't sweat over.

FEATURES

- Output rated at 55,000 LUX at 24" (610mm) distance
- 4,300° K pure white illumination
- CRI (Color Rendering Index) of 97
- Whiter, more color-correct light than halogen exam lights
- Near-zero heat from both light beam and light head housing

SPECIFICATIONS

- # of LED's: 8
- Output: 55,000 Lux @ 24"
- Focal Length: 24"
- Dimming Levels: 2
- Dimming Percentages: 100%, 50%

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)