

Luvis S200 LED Minor Surgical Light | Wall Mounted



Stock Photo - Options May Vary

FACTORY NEW

The Luvis S200 Wall Mounted LED Minor Surgical Light is designed for minor procedures and examinations in small spaces. It delivers shadow-less illumination with minimal shadow even under 50% screened conditions, while the fanless design eliminates heat and noise generation.

The S200 features an adjustable illumination pattern ranging from 18cm to 30cm, three selectable color temperature options, and intuitive controls. It delivers 100,000 lux intensity with a Color Rendering Index of 90 for accurate visualization during procedures.

Available in single or dual-head configurations, the S200 offers flexible installation options to accommodate various clinical environments and workflow requirements.

FEATURES

- Excellent Illumination
- Ergonomic & Friendly Design
- Detachable & Autoclavable Main Handle
- No UV Emission & Cool Light
- 3 Selectable Color Temperature Options
- Excellent Scalytic Effect
- Intensity: 100,000 lux
- Color Rendering Index: 90

SPECIFICATIONS

- LEDs: 18
- CRI: 90
- Color Temperature: 3,800K /4,300K/4,800K
- Head Size: 13.78" x 13.39"
- Focus Size: 7" (Circle)/11.8" (Circle)
- Intensity: Max. 100,000lux
- Power Consumption: 25W
- LED Life Span: 50,000hr

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)