

Burton AIM LED Wall Model | ALEDW



Stock Photo - Options May Vary

FACTORY NEW

The Burton AIM LED Wall Model ALEDW is a diagnostic examination light featuring state-of-the-art LED technology designed for superior illumination and energy efficiency. It delivers 45,000 Lux at 1 meter with a 92 CRI rating for accurate color visibility essential to precise diagnosis. The LED light source provides 20,000 hours of lamp life, reducing bulb replacement costs by over \$1,000 compared to traditional lighting.

The AIM LED uses 55% less energy than conventional surgical lights while maintaining exceptional performance. The removable sterilizable handle allows adjustment of the light pattern, and the unique Y-shaped design facilitates optimal positioning around the physician's head. Multiple reflectors can be tilted by simply rotating the sterilizable light handle.

The ceiling-mounted arm system offers 360° limitless rotation with a maximum reach of 63 inches and a 10-foot diameter coverage area. The light field diameter ranges from 13 to 15 inches with a depth of illumination of 54 inches. The unit is manufactured in the USA and backed by a 5-year limited warranty.

FEATURES

- Efficient LED Light Source
- On-Hand Quick Release for Removable Handle
- Removable Sterilizable Handle to Adjust Light Pattern
- Larger and Adjustable Light Pattern
- Multiple Reflectors can be tilted simply by Rotating the Sterilizable Light Handle
- Unique Y-shaped design facilitates placement of lighthouse around the Physician's head
- 360° Limitless arm-system-rotation on single Ceiling Models
- 5-year limited warranty

SPECIFICATIONS

- Illuminance: 45,000 Lux at 1 Meter
- Color Temperature: 4,300°K
- CRI (Color Rendering Index): 92
- The diameter of the Lighthouse: is 20 "
- Rated Life Of LED Lamp: 20,000 Hours
- Swivel Radius of Lamp Housing-Ceiling Mounted: 63"
- Power: 84 Watts
- Certifications / Approvals: UL60601-1 /IEC60601-1-2 / IEC60601-2-41 MDD 93/42 EEC (CE Approved) CAN/CSA C22.2601-1 M90

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)