

Steris Amsco Quantum Refurbished Surgical Lights | SQ240



Stock Photo - Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

Steris Amsco SQ240 Quantum SQ240 Refurbished Surgical Lights. This refurbished Steris Amsco SQ240 is an Excellent Choice for Any Surgical Setting. The Quantum boasts 12,000 Foot Candles of illumination, auto-re-lamping if the primary light bulb fails—note: Suspension Tubes Custom Cut for Your Required Ceiling Height.

Steris Amsco Quantum Features: Automatic Re-Lamping, Variable Spot Adjustment, and 360 Degree Mount. Includes Intensity Controller (wall mounted) and adjustable light pattern. The Amsco Quantum is a major surgical light with all of the reliability that comes with Amsco products.

Installation is available. Have Questions? Contact us today for questions or available options to fit your specific needs!

FEATURES

- Excellent Choice for any Surgical Setting
- 10,000 Ft Candles Per Head
- Automatic Re-Lamping
- Variable Spot Adjustment
- 360 Degree Mount

SPECIFICATIONS

- 10,000 Ft Candles Per Head
- 360 Degree Mount
- CRI of 92
- Lateral Illumination: 90°
- Lamp Life: 1,000 hours

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)