

Burton AIM-100 Dual Ceiling | A100DC



Stock Photo - Options May Vary

FACTORY NEW

Burton A100DC AIM-100 Dual Ceiling. Providing both excellent value and unparalleled aesthetics, Burton's AIM-100™ Dual Ceiling A100DC also boasts limitless rotation on ceiling mounts (a feature previously available only on OR-type lights).

Highest intensity illumination, the AIM-100® delivers very smooth focus and superior maneuverability. Burton AIM-100 Dual Ceiling Features: A longer arm system ensures that the light can be positioned in the most comprehensive range of positions possible, making the AIM-100® ideal for outpatient procedures in a variety of settings and practices, including trauma centers, ERs, ambulatory surgery centers, ICUs, physician offices, clinics, plastic surgery, dermatology, and ophthalmology.

FEATURES

- High intensity of 14,800 fc (160,000 lux) at 24" and 5800 fc (63,000 lux) at 1 meter
- 3500° K color temperature offers natural and true coloring
- Uses three 35-watt quartz halogen bulbs – each with a 2,000-hour rated bulb life. Bulbs are wired in parallel for continuous operation.
- Compact 20" (51 cm) diameter light head
- Extremely smooth focus mechanism with a center Sterilizable handle.

SPECIFICATIONS

- Illuminance: * 14,800 FC (160,000 Lux) at 24" and 5800 fc (63,000 Lux) at 1 meter
- Color temperature: 3500° K
- The diameter of the light head: is 20" (51 cm)
- Light field diameter: 9.6" - 12" (25 cm - 31 cm)
- The depth of Illumination: 47" (1.2 m)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)