

# Medical Illumination MI-550 Series LED Wall Mount | XLDE-WM



*Stock Photo - Options May Vary*

## FACTORY NEW

The MI-550 Series LED Wall Mount is designed for surgical procedures, clinical examinations, and diagnostic work where precise visualization is critical. The light delivers 55,000 LUX at 24 inches with 4,300°K pure white illumination and a Color Rendering Index of 97, providing superior color accuracy compared to traditional halogen exam lights. The compact LED design generates virtually no heat, operates at only 16W power consumption, and requires minimal maintenance due to long-lasting LED technology.

The MI-550 offers exceptional brightness and color fidelity with remarkable shadow control, enabling clinicians to work comfortably for extended periods. The light head weighs just 5.6 lbs and features smooth articulation with +15° upward and -60° downward vertical movement, plus approximately 540° horizontal rotation capability. Two dimming levels (100% and 50%) allow adjustment for different clinical applications and patient comfort.

This wall-mounted examination light combines high performance with energy efficiency, consuming approximately half the electricity of standard halogen lights while delivering superior illumination quality. The durable construction and efficient LED technology result in lower overall maintenance costs and extended operational life.

## FEATURES

- Output rated at 55,000 LUX at 24" (610mm) distance
- 4,300° K pure white illumination
- CRI (Color Rendering Index) of 97
- Whiter, more color-correct light than halogen exam lights
- Near-zero heat from both light beam and light head housing
- High efficiency with half the electric consumption of standard halogen lights
- Compact head design with smooth arm articulation
- Lower overall maintenance costs due to long-lasting LED's

## SPECIFICATIONS

- # of LED's: 8
- Output: 55,000 Lux @ 24"
- Focal Length: 24"
- Supply Circuit: 100-240 VAC, 50/60Hz
- Light Head Power Consumption: 16W
- Light Head Weight: 5.6 lbs
- Light Head Dimensions: 9.2" Dia x 2.6"H
- Articulating Arm Vertical Movement: +15° Up, -60° Down

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)