

Burton AIM-200™ Double Ceiling | A200DC

FACTORY NEW

Burton A200DC – AIM-200™ Double Ceiling. Philips Burton's new AIM-200 OR light provides versatility, amazing illumination, unparalleled aesthetics, and features, making it the right light for any procedure.

Additionally, with a true rendition of colors, this light is ideal for long hours of use with three cool 35-watt IRC halogen bulbs – each rated with a 1400-hour average bulb life. Burton AIM-200DC Features The right light for Ambulatory Surgical Centers.

Burton's new AIM-200 Dual Ceiling OR Light offers versatility, exceptional illumination, and unmatched aesthetics for any procedure. Moreover, it features a Color Rendering Index (CRI) of 96, ensuring accurate color representation.



Stock Photo – Options May Vary

FEATURES

- Seismic Calculations Available
- Y-shaped light head for Easy Positioning Around the Physician's Head
- Easy Bulb Replacement System From the Top of the Light head
- Removable, Sterilizable Handle to Adjust the Light Pattern
- One-Hand Quick Release for Removable Handle

SPECIFICATIONS

- Illuminance: * 90,000 Lux (8300 FC) at 1 meter (39.37 in.)
- Color temperature: 3700 k
- CRI: 96
- The diameter of the Light Head: is 20 in.
- Light Field Diameter: 10 in. - 12 in.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)