

Burton AIM-200™ Single Ceiling | A200SC



Stock Photo - Options May Vary

FACTORY NEW

Burton A200SC – AIM-200™ Single Ceiling. Philips Burton's new AIM-200 OR provides versatility, amazing illumination, unparalleled aesthetics, and features, making it the right light for any procedure.

The true rendition of colors, this light is ideal for long hours of use with three cool 35-watt IRC halogen bulbs – each rated with a 1400-hour average bulb life. Burton AIM-200 Features: The AIM-200 OR is the first Operating Room light with floor, wall, and single and double-ceiling mounts available.

With its 360° limitless arm-and-mounting system rotation around the vertical axes on the single ceiling models, the AIM-200 OR is the most versatile and practical light for the Operating Room environment.

FEATURES

- Y-shaped light head for Easy Positioning Around the Physician's Head
- Easy Bulb Replacement System From the Top of the Light head
- Removable, Sterilizable Handle to Adjust the Light Pattern
- One-Hand Quick Release for Removable Handle
- Stand-By On/Off Switch

SPECIFICATIONS

- Illuminance: * 90,000 Lux (8300 FC) at 39.37 in.
- Color temperature: 3700 k
- CRI: 96
- The diameter of the Light Head: is 20 in.
- Light Field Diameter: 10 in. - 12 in.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)