

Burton AIM-200™ Floorstand | A200FL



Stock Photo – Options May Vary

FACTORY NEW

Burton A200FL – AIM-200™ Floorstand. Philips Burton's new AIM-200 OR provides versatility, amazing illumination, unparalleled aesthetics, and many more features that make it the right light for any procedure. Delivering a true rendition of colors, this light is ideal for long hours of use with three cool 35-watt IRC halogen bulbs – each rated with a 1400-hour average bulb life.

Burton AIM-200 Floorstand Features: The AIM-200 OR is the most versatile and practical light for the Operating Room environment. Additionally, the diameter of the light head is 51 cm (20 in.). Furthermore, the AIM-200FL has a light field diameter ranging from 10 in.

to 12 in.

FEATURES

- Y-shaped light head for Easy Positioning Around the Physician's Head
- Easy Bulb Replacement System From the Top of the Light Head
- Removable, Sterilizable Handle to Adjust the Light Pattern
- One-Hand Quick Release for Removable Handle
- Stand-By On/Off Switch

SPECIFICATIONS

- Illuminance: * 90,000 Lux (8300 FC) at 39.37 in.
- Color temperature: 3700 k
- CRI: 96
- The diameter of the Light Head: is 20 in.
- Light Field Diameter: 10 in. - 12 in.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)