

Burton AIM-50 Single Ceiling Mount | A50SC



Stock Photo - Options May Vary

FACTORY NEW

Burton A50SC - AIM-50 Single Ceiling Mount. Providing both amazing value and unparalleled aesthetics, Philips Burton's AIM-50 also boasts limitless rotation on ceiling mounts (a feature previously available only on OR-type lights).

Along with the highest intensity illumination, AIM-50 delivers very smooth focus and superior maneuverability. Burton AIM-50 Single Ceiling Mount Features: A long arm system for versatile positioning, perfect for examinations.

With a CRI of 96, 45,000 Lux at 36 inches, and a 3100° Kelvin color temperature, it ensures accuracy. It uses three 50-watt quartz halogen bulbs, each with a 2000-hour average bulb life for long-term use.

FEATURES

- Illumination: 99,000 Lux at 24", 45,000 Lux at 36 inches
- Color temperature: 3100° K
- CRI (Color Rendering Index): 96
- Compact 20-inch/51cm diameter lighthead
- A larger and adjustable light pattern

SPECIFICATIONS

- High-intensity of 45,000 lux (4181 fc) at 1 meter
- High CRI (Color Rendering Index) of 96
- Uses three quartz halogen bulbs - each rated with a 2000-hour bulb life<
- Bulbs wired in parallel for continuous operation
- Compact 20" (51 cm) diameter lighthead

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)