

Burton Super Exam LED Table Mount | SELEDT



Stock Photo - Options May Vary

FACTORY NEW

The Burton SELEDT Super Exam LED Table Mount delivers powerful, focused illumination for medical exams and procedures. It produces up to 35,000 lux with a 4000K color temperature and CRI above 95 for accurate tissue and skin tone rendering. This ensures consistent lighting performance across exam rooms and specialty environments.

Designed for flexibility, this table-mounted exam light clamps securely to flat surfaces up to three inches thick. The drift-free K-arm allows smooth positioning without unwanted movement, and users can easily adjust the light with one hand during procedures. An ergonomic handle integrates power and dimming controls for quick access and improved workflow efficiency.

Built for durability, the Super Exam LED features a 50,000-hour LED lifespan and energy-efficient 20-watt operation. The SteriTouch® coating helps inhibit bacterial growth for improved hygiene, and the smooth exterior supports fast, thorough cleaning between patients. This delivers reliable performance and long-term clinical value.

FEATURES

- High-performing energy-efficient LED module with 50,000-hour rated life
- Flexible drift-free spring-balanced arm
- Power and dimming switch positioned near the handle
- SteriTouch® coating inhibits the growth of bacteria
- 35,000 lux (3251 fc) at 39.37" (1 m)
- 4,000 K color temperature
- CRI (Color Rendering Index) of >95
- Industry-leading 5-year warranty

SPECIFICATIONS

- Illuminance: 35,000 lux (3251 fc) at 39.37" (1 m)
- Color temperature: 4000 K
- CRI >95
- Light sources: Light Emitting Diode, minimum 50,000-hour rated life
- Power: 20 Watts
- Mounting options: Single Ceiling, Table, Floor and Wall
- Applications: General examination
- Product warranty: 5-year parts

**Brochure and images are a general overview of the product. Actual product options may vary.*