

# Clinton 48"L Maple Therapy Table w/Hand Crank Height Adj | 7445M

## FACTORY NEW

Clinton 48"L Maple Therapy Table w/ adj Height, Model 7445M (previously 77-45CM). Clinton Industries is a true American success story of entrepreneurial fortitude and industry innovation. Started in 1979 by founder and present owner, Robert C.



Bohn, with a singular focus to build a better more affordable treatment table, Clinton Industries flourished from a "garage start-up" with a handful of treatment tables to a multi-million dollar industry. The Clinton 7445M Therapy Table features a classic maple butcher block top, easy height adjustment with a fold-away handle, an aluminum and powder-coated base, and dependable hydraulic height adjustment.

The sturdy construction ensures durability for frequent use.

*Stock Photo - Options May Vary*

## FEATURES

- Classic maple butcher block top
- Easy height adjustment with fold-away handle
- Aluminum and powder-coated base
- Dependable, hydraulic height adjustment

## SPECIFICATIONS

- Model Number: 7444M
- Length: 48"
- Width: 36"
- Adjustable Height: 27" to 39"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free      Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)