

Clinton Maple 72"L Therapy Table w/Hand Crank Height Adj | 7444M

FACTORY NEW

Clinton Maple 72"L Therapy Table w/ Adj Height: 7444M (previously 77-44CM). Clinton Industries is a true American success story of entrepreneurial fortitude and industry innovation. Started in 1979 by founder and present owner, Robert C.



Stock Photo - Options May Vary

Bohn, with a singular focus to build a better more affordable treatment table, Clinton industries flourished from a "garage start-up" with a handful of treatment tables to a multi-million dollar industry. The Clinton 7444M Therapy Table features a classic maple butcher block top, easy height adjustment with a fold-away handle, an aluminum and powder-coated base, and dependable hydraulic height adjustment.

Its durable construction suits demanding therapy environments. The adjustable height offers versatility for various treatments.

FEATURES

- Classic maple butcher block top
- Easy height adjustment with fold-away handle
- Aluminum and powder-coated base
- Dependable, hydraulic height adjustment

SPECIFICATIONS

- Model Number: 7444M
- Length: 72"
- Width: 36"
- Adjustable Height: 27" to 39"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)