

Conmed 2450 Electrosurgical Unit - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Conmed System 2450 Electrosurgical Unit is an advanced surgical platform designed for safe, precise tissue cutting and coagulation. It delivers multiple cutting modes—Pure Cut for clean incisions and Blend for hemostasis during cutting—along with four coagulation options for effective dissection and bleeding control. Bipolar functionality includes both Micro and Macro settings for precise desiccation with minimal lateral tissue damage.

Dynamic Response Technology continuously adjusts energy output based on tissue type, reducing unwanted thermal spread. The Automatic Return Monitor (ARM) detects issues with patient grounding pad placement, adding an essential safety layer. The intuitive operator interface enables quick setup and operation in busy surgical environments.

The unit includes a monopolar foot switch and Auto Programming to store past settings. Multiple accessory ports accommodate additional equipment, with optional bipolar foot switch and rolling cart available. This reliable system provides consistent electrosurgical performance for clinics and hospitals.

FEATURES

- Smoother Cut through all Tissue Types
- 3 Cutting Modes: Low Cut, Pure Cut, and Blend for Cutting with Hemostasis
- Instant Response Technology
- Power Efficiency Rating (PER) 98
- Accurate and Consistent Cut Performance
- 4 Coag Modes
- Monopolar Foot-Switch Included
- Completely Refurbished

SPECIFICATIONS

- Monopolar Cut (Max): 300 Watts
- Monopolar Coag: 120 Watts
- Bipolar: 70 Watts
- Blend Cut: 200 Watts
- Power Requirements: 120V
- Weight: 32 Lbs (unit only)
- Dims: 14.4"W x 21.5"D x 7"H
- 4 Coag Modes

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)