

Conmed 5000 Electrosurgical Unit - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Conmed System 5000 Electrosurgical Unit is a versatile, technologically advanced surgical platform offering multiple operation modes tailored to specific procedures and surgeon preferences. Specialty modes include General, Fluids (which boosts energy in fluid environments), Laparoscopic (limiting power to reduce inadvertent burns), and Pulse Cut/Coagulation options.

The unit features multiple cutting modes—Pure Cut for scalpel cuts, Blend (1, 2, 3) for hemostasis while cutting, and Pulse Cut—along with coagulation modes (Pinpoint, Standard, Spray, Pulse) for dissection and hemostasis. Bipolar settings include both Micro and Macro modes for enhanced procedural flexibility.

Dynamic Response Technology reacts to tissue changes to minimize tissue damage. The Automatic Return Monitor (ARM) detects faults in the patient grounding circuit to help prevent pad site burns. The unit includes programmable memory to recall up to 9 user settings and comes with a monopolar foot switch.

FEATURES

- Specialty Modes: General, Fluids, Laparoscopic, Pulse Cut/Coag
- Cutting Modes: Pure, Blend(3), Pulse
- Coag Modes: Spray, Pinpoint, Standard, Pulse
- Bipolar Modes: Micro and Macro Bipolar
- Dynamic Response Technology
- Programmable Memory Settings (9)
- Monopolar Foot-Switch Included
- Completely Refurbished

SPECIFICATIONS

- Monopolar Cut (Max): 300 Watts
- Monopolar Coag: 120 Watts
- Blend Cut (1-3): 200 Watts
- Spray Coag: 80 Watts
- Micro BiPolar: 50 Watts
- Macro Bipolar: 90 Watts
- Power Requirements: 120V
- Dims: 13.5"W x 21.5"D x 5.5"H

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)