

Conmed BiCap III Electrosurgical Unit – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

The Conmed BiCap III Electrosurgical Unit delivers precision and control for advanced surgical procedures with Dynamic Response Technology that adjusts to tissue impedance in under one millisecond, ensuring consistent power and reduced tissue damage. Front panel volume control and Pulse Cut Mode enable accurate energy bursts during critical dissections, while automatic programming restores previous settings instantly for enhanced efficiency.

The BiCap Superconductor Probes feature three electrode pairs with a strong ceramic tip and large electrode area that spreads energy evenly for effective bipolar hemostasis while reducing tissue sticking. The short 3.3mm tip facilitates easier use with deflected scopes, and color-coded shafts provide clear size identification.

A central irrigation channel improves visibility during surgery for cleaner, more controlled work. The compact design and reliable performance deliver dependable results in demanding surgical environments.

FEATURES

- Pure Cut: 200 Watts
- Blend Cut: 200 Watts
- Coag: 80 Watts
- Bipolar: 50 Watts
- Front panel volume control
- Pulse Cut Mode provides precise energy bursts for critical dissections
- Automatic Programming option restores the BiCap's last settings used
- Less than 1 millisecond reaction time

SPECIFICATIONS

- Weight: 14.65 lb (6.65 kg) (without cart)
- Dimensions: 6" x 13" x 18" (15.4 cm x 33 cm x 45.7 cm)
- Modes: Pure Cut, Blend Cut, Pulse Cut, Coag, and Bipolar
- Output Control System: Dynamic response
- Contact Quality Monitoring System: ARM
- Duty Cycle: 15 seconds on, 30 seconds off
- Technology: Isolated RF Output
- Power: 120 V

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)