

Covidien Valleylab FT10 Energy Platform - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Valleylab FT10 Energy Platform represents nearly 50 years of innovation in energy-based surgical devices. It delivers enhanced performance, greater precision, and improved efficiency across a full range of surgical procedures, enabling surgeons to operate with increased speed, accuracy, and control in demanding surgical environments.

Key performance improvements include LigaSure vessel sealing up to 50% faster with lower jaw temperatures, quicker auto bipolar technology activation, and stronger monopolar performance. TissueFect sensing technology adjusts energy in real time, applying the precise amount needed exactly where and when required.

These innovations streamline surgical workflows, reduce overall procedure time, and support better patient outcomes. The FT10 Energy Platform provides reliable performance and consistent results for a wide range of clinical applications.

FEATURES

- 7 inch LCD Touchscreen
- Improved Monopolar and Bipolar Performance
- Up to 50% Faster Sealing times
- TissueFect™ Tissue Sensing Technology reads tissue composition in real time
- Closed-loop control at 434,000 data points per second (100x more than Force Triad)

SPECIFICATIONS

- Pure Mode: 1- 300 W / 1287 V
- Blend Mode: 1 -200 W / 2178 V
- Valleylab Mode: 5 - 85 W / 2783 V
- Soft Mode: 1 - 120 W / 3449 V
- Spray Mode: 1 - 120 W / 3933 V
- 14.5"W x 18.2"D x 7"H
- 22.3 lbs / 10.1 kgs

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)