

# ValleyLab Force FX-c Electrosurgical Unit – Refurbished



Stock Photo – Options May Vary

## FULLY REFURBISHED • CLINICALLY CERTIFIED

The ValleyLab Force FX-c is a fully refurbished electrosurgical unit designed for general surgery applications. It delivers reliable performance across all tissue types with smooth cutting action and precise control. The unit features durable design and trusted technology suitable for hospitals, clinics, and surgical centers.

The Force FX-c offers three cutting modes: Low Cut, Pure Cut, and Blend for cutting with hemostasis. It includes Instant Response Technology and achieves a Power Efficiency Rating (PER) of 98 for accurate and consistent cut performance. The unit also provides four coagulation modes for versatile surgical applications.

A monopolar foot switch is included with the unit. Optional bipolar foot-switch and rolling cart are available for additional functionality and mobility.

## FEATURES

- Smoother Cut through all Tissue Types
- 3 Cutting Modes: Low Cut, Pure Cut, and Blend for Cutting with Hemostasis
- Instant Response Technology
- Power Efficiency Rating (PER) 98
- Accurate and Consistent Cut Performance
- 4 Coag Modes
- Monopolar Foot-Switch Included
- Completely Refurbished

## SPECIFICATIONS

- Monopolar Cut: 300
- Monopolar Coag: 120
- Bipolar: 70

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)