

Aaron Bovie Specialist PRO Electrosurgical Generator | A1250S



Stock Photo - Options May Vary

FACTORY NEW

The Aaron Bovie Specialist PRO Electrosurgical Generator (A1250S) is a multi-purpose electrosurgical unit designed for physician offices and specialty facilities. It features a modern push-button interface capable of performing most monopolar and bipolar surgical procedures with five energy outputs delivering up to 120 watts of power across Cut, Blend, Coagulation, Fulguration, and Standard Bipolar modalities.

The unit incorporates advanced safety features including BovieNEM™ (Neutral Electrode Monitoring) pad-sensing technology that monitors the return electrode for optimal patient protection, BovieDED™ (Digital Error Detection) that monitors all generator functions, audible alert tones, and self-test circuits. The design provides large illuminated digital displays and ample raised push-button power controls for rapid adjustments.

The Specialist PRO operates with the reliability and functionality of an operating room generator in a compact, lower-powered package suitable for smaller facilities. It is backed by Bovie's four-year manufacturer's limited warranty and includes isolated circuitry and discrete output management.

FEATURES

- Designed to operate like an OR generator in a smaller, lower-powered package
- State-of-the-art safety features including Bovie NEM (Neutral Electrode Monitoring) pad-sensing technology, which monitors the return electrode for optimum patient protection
- BovieDED (Digital Error Detection) monitors all functions of the generator to ensure constant safety for the patient, doctor, and staff
- Backed by Bovie's four (4) year manufacturer's limited warranty
- Large illuminated digital displays
- Ample raised push button power controls provide rapid power adjustments
- Self-test circuits
- Audible tones

SPECIFICATIONS

- Cut (120 Watts)
- Blend (90 Watts)
- Coagulation (80 Watts)
- Fulguration (40 Watts)
- Bipolar (30 Watts)
- Monopolar Footswitch
- Monopolar Hand Switching Receptacle Accepts standard three-pin handpiece
- Bipolar Receptacle Accepts standard cable for bipolar handpieces

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com