

# Aaron/Bovie High-Frequency Electrosurgery Generator | 1250



Stock Photo - Options May Vary

## FACTORY NEW

The Aaron 1250 by Bovie is a multipurpose electrosurgical generator designed for use in physician offices and surgical centers. It provides monopolar and bipolar functionality with five output modes: Cut (120 watts), Blend (90 watts), Coagulation with pinpoint and fulguration settings (80 and 40 watts respectively), and Bipolar (30 watts). The unit delivers flexible surgical capabilities with reliable performance.

Safety features include self-test circuits, audible tones, isolated circuitry, and Bovie NEM (Neutral Electrode Monitoring) for automatic protection during procedures. The generator includes a large illuminated digital display and a rotary power-control dial for smooth, rapid adjustments during surgery.

The unit comes complete with a power cord, sterile disposable pencil, reusable grounding cord, sterile disposable split grounding pads, and three sterile disposable electrodes in various configurations. A four-year warranty is included with purchase.

## FEATURES

- Large illuminated digital displays
- The rotary power-control dial provides smooth, rapid power adjustments
- Five Output Modes: Cut (120 Watts), Blend (90 Watts), Coagulation (Pinpoint 80 Watts), Fulguration (Spray 40 Watts), Bipolar (30 Watts)
- Self-test circuits
- Audible tones Discreet outputs (FCFS™ - First Come First Served)
- Isolated circuitry BovieNEM™ (Neutral Electrode Monitoring)
- Digital Error Detection
- Included: Unit, Power Cord, One Sterile Disposable Pencil (ESP1), One Reusable Grounding Cord (1252C), Five Sterile Disposable Split Grounding Pads (ESRE-1), Three Sterile Disposable Electrodes (one ES01 blade, one ES02 needle, one ES20 ball), and User's Guide

## SPECIFICATIONS

- Output Power Cut: 120 watts
- Blend: 90 watts
- Coagulation: 80 watts
- Fulguration: 40 watts
- Bipolar: 30 watts
- Line Voltage A1250: 120 VAC ± 10% / A1250-220: 240 VAC ± 10%
- Output Frequency 310-550 kHz
- Dimensions and Weight Width: 26.04 cm (10.25") Height: 15.24 cm (6") Depth: 29.85 cm (11.75") Weight: < 4 kg (< 9 lbs.)

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)