

Covidien ValleyLab Force Triad Electrosurgical Unit – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Covidien ValleyLab Force Triad is a comprehensive electrosurgical unit that combines monopolar, bipolar, and vessel sealing capabilities in a single system. It delivers precisely controlled energy for improved tissue effects across multiple surgical applications.

The unit features a user-friendly interface with three touchscreen LCD displays that simplify navigation and menu selection. Backward compatibility with common electrosurgical accessories allows integration with existing operating room inventory.

The ValleyLab mode provides monopolar hemostasis with dissection at lower power settings, resulting in reduced charring, thermal spread, and arcing for smoother tissue passage. The system includes monopolar and LigaSure footswitches and incorporates TissueFect sensing technology and REM ground pad contact monitoring.

FEATURES

- Full Features of an ESU with the LigaSure Vessel Sealing System
- 5 Modes: Cut, Valleylab, Coag, Bipolar, and LigaSure
- 3 LCD Screens, all Touchscreens
- TissueFect™ Sensing Technology
- REM Ground Pad Contact Monitoring
- 4 Coag Modes
- Monopolar and Ligasure Footswitches Included
- Completely Refurbished

SPECIFICATIONS

- Cut: 660 Hz
- ValleyLab Mode: 800 Hz
- Coag: 940 Hz
- Bipolar: 940 Hz
- LigaSure: 440 Hz
- 18"W x 20"D x 10"H
- 30 Lbs

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)