

ValleyLab Force 2 Electrosurgical Unit – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

The Refurbished ValleyLab Force 2 Electrosurgical Unit is a multi-use ESU recognized worldwide for precision and reliability. It delivers powerful performance with 300W monopolar cut, 120W coagulation, and 70W bipolar capabilities, along with three preset blend modes for versatile surgical applications.

The unit incorporates enhanced safety features including REM Return Electrode Monitoring to help prevent electrode site burns. An easy-to-read digital display shows all settings, and audible alarms alert users to any operational issues.

This completely refurbished unit comes with a monopolar foot switch, power cord, and manual. Optional accessories including a cart and bipolar foot switch are available separately.

FEATURES

- Multi-Use, Powerful Unit
- 300W Cut / 120W Coag/ 70W Bipolar
- Enhanced Safety with REM Return Electrode Monitoring
- Audible Alarms
- 3 Preset Blend Modes
- Digital Display
- Includes Monopolar Foot-switch
- Refurbished

SPECIFICATIONS

- Monopolar Cut: 300 Watts
- Monopolar Coag: 120 Watts
- Bipolar: 70 Watts
- Blend 1: 250 Watts
- Blend 2: 200 Watts
- Blend 3: 150 Watts
- 21"D x 13"W x 8"H
- Weight: 23 Lbs (Unit Only)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)