

Valleylab LigaSure Vessel Sealing Unit – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

The LigaSure vessel sealing system provides permanent vessel occlusion through revolutionary technology that fuses vessel walls to create a permanent seal. This approach replaces multiple hemostatic tools by delivering reliable vessel sealing without relying on traditional ligation methods.

The unit features a digital display, bipolar and macrobipolar operating modes, and delivers constant, consistent output with an average seal cycle of 2-4 seconds. The system operates at 473 kHz across all modes for precise hemostatic control.

This refurbished unit includes a footswitch and power cord, with dimensions of 40.6 cm length × 38.6 cm width × 12.7 cm height and weighs 5.9 kg. It requires standard 120V power and comes with a 6-month warranty.

FEATURES

- Digital Display
- Bipolar Mode
- Macrobipolar Mode
- Constant, Consistent Output
- Average Seal Cycle is 2-4 Seconds
- Foot-Switch Included
- Completely Refurbished

SPECIFICATIONS

- Macrobipolar: 473 kHz sinusoid
- Bipolar: 473 kHz Sinusoid
- Seal Mode: 473 kHz Sinusoid, Pulsed
- Height (unit only): 12.7 cm (5")
- Width (unit only): 38.6 cm (15.2")
- Length (unit only): 40.6 cm (16")
- Weight (unit only): 5.9 kg (13 lbs)
- Power Requirements: 120V

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)