

Aaron Bovie 950 High-Frequency Electrosurgical Unit | A950



Stock Photo - Options May Vary

FACTORY NEW

The Aaron 950 is a dual-function electrosurgical unit combining a high-frequency desiccator with a Bovie generator. The Bovie generator provides up to 60 watts of cutting power with three modes: Cut, Blended Cut, and Pinpoint Coagulation, suitable for office-based procedures including LLETZ/LEEP. The desiccator side delivers 35 watts for fulguration and 30 watts for bipolar work.

The unit features user-friendly controls with color-coded membrane pads, a large rotary dial for mode selection, and digital error detection that continuously monitors output and disables the unit if issues arise. Six preset buttons allow one-touch recall of favorite settings. The ergonomic design includes a compact front panel and isolated power configuration.

The Aaron 950 comes with an autoclavable 3-button handpiece that rests safely in a holder when not in use. All industry-standard 3/32" electrodes are compatible. The nonabsorbent ABS enclosure and liquid-resistant membrane switches facilitate quick post-procedure cleanup. The desiccator operates without requiring a return pad.

FEATURES

- Bovie® Generator
- 60 Watts - Cut, Blend, Coagulation
- Isolated Power
- 6 Presets - For The Doctor's Favorite Settings and one-touch recall
- High-Frequency Desiccator: 35 Watts - Fulguration, 30 Watts - Bipolar, 1-10 Watts - Adjustable in 1/10 watt increments
- The desiccator does not require the use of a return pad
- Change Power In One Step, Not Four
- The Bovie 3-Button Hand-piece - it's autoclavable, so you can have sterility when you need it

SPECIFICATIONS

- Generator: Cut: 60 watts, Blend: 60 watts, Coagulation: 60 watts
- High-Frequency Desiccator: Fulguration: 35 watts, Bipolar: 30 watts
- Line Voltage: A950™: 120 VAC ± 10%
A950™-220: 230-240 VAC
- Output Frequency: 350-800 kHz
- Line Frequency: 50-60 Hz
- Dimensions: 10.25"W x 7"H x 4.5"D
- Weight: < 8 lbs.
- Warranty (Unit): 2 years

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com