

Bovie DERM 942 High-Frequency Desiccator | A942



Stock Photo - Options May Vary

FACTORY NEW

The Bovie Derm 942 High Frequency Desiccator delivers precise control for dermatologic and surgical procedures with 40 watts of monopolar and bipolar power, adjustable in 1/10th watt increments up to 10 watts. This precision allows medical professionals to desiccate and fulgurate tissue with accuracy and efficiency.

The unit features a 3-button power control handpiece or rotary dial for effortless power adjustment during use, with a large, bright LED display for excellent visibility from anywhere in the procedure room. An optional footswitch can be connected to enable foot control activation, eliminating the need for a separate handpiece and improving ergonomics.

Engineered for durability and reliability, the Derm 942 is designed for long-term performance in clinics and hospitals. It offers unmatched efficiency and control for both minor dermatologic treatments and advanced electrosurgical tasks.

FEATURES

- 40 watts of power for high-frequency desiccation
- Fine control: Adjust in 0.1-watt steps (0-10 W range in bipolar mode)
- Global use: Works with 100-240 VAC power
- No grounding pad needed (ground reference unit)
- 3-button handpiece: (1) Works with or without the optional footswitch (2) Lets you adjust power and activate directly from the pencil
- Plug-and-play footswitch (no handpiece change required)
- Built-in safety: Automatically shuts off and shows error codes if an issue occurs
- Reliable A902 handpiece: (1) 2-pin design with gold-plated collet for strong conductivity
Comfortable, (2) non-slip grip

SPECIFICATIONS

- Max Wattage: 40
- Bipolar Mode: Yes
- Voltage: 110-240V
- Electrodes Size: 3/32"
- Dimensions: 7.4"H x 8.98"W x 4.13"D
- Weight: 5 lbs
- Warranty: 4-year manufacturer's limited warranty

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com