

Bovie® Bantam PRO™ 50W Electrosurgical Generator Gynecology Package | A952-G



Stock Photo - Options May Vary

FACTORY NEW

The Bovie A952-G Bantam PRO Electrosurgical Generator combines a High Frequency Desiccator with a Bovie Generator, delivering 50 watts of cutting and coagulation power for a wide range of clinical and surgical applications. The unit features three operating modes—Cut, Blended Cut, and Pinpoint Coagulation—with easy-touch, color-coded controls and an ergonomic handpiece for procedural precision and safety.

Engineered for reliability, the A952-G uses advanced chip technology to reduce internal heat and extend equipment life. Its non-absorbing, high-impact ABS enclosure simplifies cleaning and maintenance while the compact, mobile design fits seamlessly into any clinical workspace.

The A952-G PRO Gynecology Package is specifically designed for gynecology practices and includes GYN electrodes, a mobile stand, a smoke evacuator, and essential supplies. Featuring Bovie NEM™ pad sensing technology, the unit enhances patient protection during electrosurgical procedures with continuous monitoring and digital error detection.

FEATURES

- Use of the Bovie Bantam PRO, A952, promotes procedural confidence. The handpiece fits comfortably in your hand
- When an electrosurgical generator mode, the unit features Bovie NEM™ pad-sensing technology, which monitors the return electrode for optimum patient protection
- Post-procedure cleanup is simple (and fast) due to nonabsorbent ABS enclosures and our use of liquid-resistant membrane switches that can be wiped clean with a damp cloth
- High-Frequency Desiccator 40 Watts - Fulguration 40 Bipolar 40 Micropolar .1-10 Watts - Adjustable in 1/10 watt increments
- Change Power in one step, not four - the A952 allows the doctor to change the power setting with one action, using the up/down buttons on the autoclavable handpiece or the rotary power-control dial. This saves time and protects the sterile field during sensitive procedures
- The Bovie 3-Button Handpiece is better by design - the A902 handpiece (Patent #8,998,899) allows you to adjust the power output up and down as well as to activate the unit, even when you are in foot control mode. Designed with a 2-pin plug design, no more flying lead
- Digital error detection means unsurpassed safety for you and your patient. The Bovie Bantam | PRO, A952, continuously monitors every aspect of the unit's output

SPECIFICATIONS

- Output Power: Cut-50 watts, Blend: 50-watts, Coagulation: 50-watts, Fulguration: 40-watts, Bipolar: 40-watts
- Output Frequency: 304-451 kHz
- Line Voltage: 100-240 VAC
- Line Frequency: 50-60 Hz
- Dimensions and Weight Width: 27.9 cm (11") Height: 15.2 cm (6") Depth: 39.4 cm (15.5")
- Weight: 4.3 kg (10 lbs.)
- Contact Monitoring Systems: BovieNEM™
- Cooling System: Passive Cooling

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com