

Bovie® Bantam PRO™ 50W Electrosurgical Generator | A952



Stock Photo - Options May Vary

FACTORY NEW

The Bovie A952 Bantam PRO Electrosurgical Generator combines a 50-watt electrosurgical generator with high-frequency desiccation capability in a compact unit designed for office-based procedures. Advanced chip technology and reduced internal components minimize heat for long-term reliability. The high-impact ABS enclosure resists stains and absorbs no fluids, simplifying cleanup after procedures.

The A952 offers three operating modes—Cut, Blended Cut, and Pinpoint Coagulation—providing precise control for desiccation, cutting, and coagulation. Color-coded membrane pads, a large rotary dial, and compact front panel enable seamless operation with rapid power adjustments. Settings can be fine-tuned directly from the autoclavable handpiece for added convenience.

The ergonomic handpiece includes a secure holder to minimize contact with the unit and maintain sterile field integrity. Bovie NEM™ pad sensing technology monitors the return electrode during use for patient protection. The system supports all standard 3/32" electrodes and offers optional foot control without requiring a separate handpiece.

FEATURES

- Cut Mode: Smooth, consistent tissue cutting with minimal bleeding
- Blend Mode: 55/45 mix of cutting and coagulation for improved hemostasis
- Coag Mode: Precise pinpoint coagulation for controlled bleeding
- Dual-function system with 50 watts of cutting/coagulation power and high-frequency desiccation
- Simple, color-coded controls and large rotary dial for easy adjustments
- Compatible with all 3/32" disposable or reusable electrodes; handpiece is autoclavable
- Bovie NEM™ pad sensing technology for optimal patient protection
- Optional footswitch, plug-and-play, no special handpiece needed

SPECIFICATIONS

- Output Configuration: 304-451 kHz
- Height: 7.40"
- Width: 8.98"
- Depth: 4.13"
- Weight: 5 lbs
- Input Characteristics: 100 - 240 VAC
- Output - Cut Mode: 50 W @ 500 O
- Output - Coagulation Mode: 50 W @ 1000 O

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com