

Aaron Bovie High-Frequency Desiccator Unit w/Power Control Handpiece | A940



Stock Photo - Options May Vary

FACTORY NEW

The Aaron A940™ is a high-frequency desiccator unit designed for maximum reliability and durability in electrosurgical applications. Advanced chip technology and reduced internal components minimize heat generation, extending the unit's operational life. The non-absorbing, non-staining high-impact plastic case simplifies cleaning and maintenance.

The A940™ features digital error detection that continuously monitors all output parameters for enhanced safety. If any anomaly is detected, the unit instantly disables output and displays an error code. The system delivers up to 40 watts of power with precise control in 1/10th-watt increments in both monopolar and bipolar modes.

Designed for user efficiency, the A940™ reduces procedural steps and minimizes required accessories. Universal power supply compatibility (100-240 VAC) makes it suitable for global use. The unit includes a power control handpiece, disposable dermal tips, wall mount kit, and comprehensive user documentation.

FEATURES

- More Power with 40 watts
- Higher Resolution in Bi-Polar Mode increase and a decrease in 1/10th-watt increments from zero to ten watts
- Universal Power with 100-240 VAC...the A940 is a global option
- More User-Friendly with features that can save you time and money
- Less to Do - fewer Steps to Accomplish the Same Procedures
- Less to buy, Fewer accessories needed for additional functionality
- A940 Unit, Power Cord
- 2-Year Warranty

SPECIFICATIONS

- Height: 15.9
- Width: 22.2
- Depth: 11.4
- Weight: < 3 lbs.
- Output Power: To 40 watts (0-10 watts in 1/10th-watt increments on both mono and bipolar)
- Line Voltage: Universal Power Supply A940: 100 to 240 VAC
- Output Frequency: 550 kHz
- Line Frequency: 50 - 60 Hz

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)