

Aaron Smoke Shark II™ w/Long Life 35 Hour Filter & Pneumatic Footswitch | SEO2



Stock Photo - Options May Vary

FACTORY NEW

The Smoke Shark II is a modern, lightweight smoke evacuator with a compact design that simplifies use, storage, and stackability. It comes packaged with the SF35 Extended-Life Filter and FSSEP Pneumatic Footswitch, eliminating the need for separate orders.

The SF35 Extended-Life Filter provides up to 35 hours of performance, reducing procedure costs by approximately 38% compared to the previous model when used at the lowest setting. The covered 3-port filter design accommodates tubing sizes of 1/4", 3/8", and 7/8", providing versatility for various procedure setups and accessories.

The unit delivers excellent smoke evacuation performance with an ultra-quiet design. It is -3 and RoHS compliant, meeting regulatory standards for medical equipment use.

FEATURES

- Easy To Operate
- Quiet And Effective Smoke Evacuation
- Long-life 35-hour Filter
- Low Cost Per Procedure
- Excellent Particle Capture
- Includes: Unit, Power Cord, User's Guide. Pneumatic Footswitch (FSSEP) and Long Life 35 Hour Filter (SF35)

SPECIFICATIONS

- Line Voltage SE02: 100/120 VAC, SE02-220: 220/240 VAC
- Line Frequency 50-60 Hz
- Width: 11"W x 6"H x 15.5"D
- Weight: 10 lbs.
- Warranty (Unit) - 2 year

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)