

# Bovie Derma 101 High-Frequency Desiccators | D101



Stock Photo - Options May Vary

## FACTORY NEW

The Bovie Derma 101 is a 10-watt high-frequency desiccator designed for in-office minor skin procedures including removal of skin tags, benign lesions, and premalignant lesions. The unit delivers monopolar power adjustable in 1/10th-watt increments for precise tissue desiccation and fulguration.

Power adjustment is available via either a 3-button handpiece control or rotary dial on the unit. A large, bright LED display provides excellent visibility during procedures. Optional foot control compatibility is available with the A803 Bovie Footswitch (sold separately).

The unit arrives packaged and ready to use, complete with 18 disposable dermal tips. The economical design allows practitioners to recover the equipment cost through procedure revenue while the included tips are consumed.

## FEATURES

- 10-watt high frequency
- Ideal for performing in-office minor skin procedures such as the removal of skin tags, benign lesions, and premalignant lesions
- Monopolar power: Adjustable in 1/10th-watt increments to desiccate and fulgurate tissue
- Adjustable Power: Adjusted using the 3-Button Power Control Handpiece (A902) or using the rotary dial on the unit
- Unit Display: Large, bright, LED display for excellent visibility across the procedure room
- Optional A803 Bovie Footswitch (sold separately) compatibility
- The unit is packaged and ready to use
- Comes with a total of 18 disposable dermal tips

## SPECIFICATIONS

- Height 2.7" (69mm) Width 7.2" (184mm) Length 8.2" (208mm)
- Weight < 1.4 kg (< 3 lbs)
- Main Voltage: 100 - 240 VAC  $\pm$  10%
- Main Frequency: 50 - 60 Hertz Main Current: 320mA Max
- Power Consumption: 28 VA
- Mode of Operation: Intermittent Operation
- Duty Cycle: 10sec on / 30sec off
- four years Four-Year Warranty

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)